

Product Datasheet

Canine Oncostatin M Recombinant Protein (orb1820210)

Description	The Canine Oncostatin M (OSM) yeast-derived recombinant protein is not tagged and is naturally endot...
Target	Oncostatin M
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	SCSDKYPELL GQLQKQADFM QHTNTLLDLY IRSQGLDKNG LKEHCRERPG AFPSKDALQR LSRRVFLRTL DTTLGQVLLR LAALEQDIPK AQDLEMLSGV KLNIRGFKNN IHCMAQLLPG SSETTEPTPT SPGASPSPTP TLDTFQRRLE GCRFLHGYHR FMRSVGQVFR EWGKSLRSR R (181)
Purity	98%
Source	Yeast
MW	20.6 kDa
Entrez	611921
Expiration Date	6 months from date of receipt.