

## Product Datasheet

### Rabbit IL-21 Recombinant Protein (orb1820220)

<b>Description</b>	The Rabbit IL-21 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. Th...
<b>Target</b>	IL-21
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	HKSSSKGQDR YMIRMHQLLD IVDQLQSDVN DLDPDFLPAP QDVQKGCEQS AFSCFQKAQL KPANAGDNGK RISSLIKQLK RKL PSTKSKK TQKHRPTCPS CYSYEKKNLK EFLERLKSLI QKMIHQHLL EHLR (133)
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	15.4 kDa
<b>Entrez</b>	<a href="#">100344172</a>
<b>Expiration Date</b>	6 months from date of receipt.