

Product Datasheet

Rabbit IL-17F Recombinant Protein (orb1820229)

Description	The Rabbit IL-17F yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. T...
Target	IL-17F
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	ARKNPKVEDA TLLKPGNCPD LEDNSVRVDI RILSQNHGVP ITHDFQNRSI SPWDYNITRD PHRFPSEIAE ARCRHSGCIN AQGEEDSSMN SVPIQQEFLV LQREPQGCSR SFRLEKVVWT VGCTCVTPIV RYAA (134) [Arg102Gln]
Purity	98%
Source	Yeast
MW	15.1 kDa
Entrez	100339570
Expiration Date	6 months from date of receipt.