

Product Datasheet

Human IL-17F Recombinant Protein (orb1820232)

Description	The Human IL-17F yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. Th...
Target	IL-17F
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	ARKIPKVGHT FFQKPESCPP VPGGSMKLDI GIINENQRVS MSRNIESRST SPWNYTVTWD PNRYPSEVVQ AQCRNLGCIN AQQKEDISMN SVPIQQETLV VRRKHQGCSV SFQLEKVLVT VGCTCVTPVI HHVQ (134)
Purity	98%
Source	Yeast
MW	15.0 kDa
Entrez	112744
Expiration Date	6 months from date of receipt.