

Product Datasheet

Feline IL-5 Recombinant Protein (orb1820234)

Description	The Feline IL-5 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
Target	IL-5
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	IAVQSPMNRL VAETLALLST HRTLLIGDGN LMIPTPEHNN HQLCIEEVFQ GIDTLKNRTV PGDAVEKLFNR NLSLIKEHID RQKKKCGGER WRVKKFLDYL QVFLGVINTE WTMS (115)
Purity	98%
Source	Yeast
MW	13.2 kDa
Entrez	493803
Expiration Date	6 months from date of receipt.