

## Product Datasheet

# Ovine IL-12 p40 Recombinant Protein (orb1820236)

<b>Description</b>	The Ovine/Caprine IL-12/IL-23 p40 yeast-derived recombinant protein is not tagged and is naturally e...
<b>Target</b>	IL-12 p40
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	IWELEKNVYV VELDWYPNAP GETVVLTCDT PEEDGITWTS DQSSEVLGSG KTLTIQVKEF GDAGQYTCHK GGEVLSRSL LLHKKEDGIW STDILKDQKE PKAKSFLKCE AKDYSGHFTC SWLTAISTNL KFSVKSSRGS SDPRGVTCGA ASLSAEKVSM DHREYNKYTV ECQEGSACPA AEESLPIEVV MEAVHKLKYE NYTSSFFIRD IIKPDPPKNL QLRPLKNSRQ VEVSWEYPDT WSTPHSYFSL TFCVQVQGKN KREKKLFTDQ TSAKVTCHKD ANIRVQARDR YYSSFWSEWA SVSCS (305)
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	34.5 kDa
<b>Entrez</b>	<a href="#">443472</a>
<b>Expiration Date</b>	6 months from date of receipt.