

Product Datasheet

Swine C3a Recombinant Protein (orb1820245)

Description	The Swine C3a yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The p...
Target	C3a
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	SVQLMEKRMD KLGQYSKDVR RCCEHGMRDN PMKFSCQRRRA QFIQHGDACV KAFLDCCEYI AKLRQQHSRN KPLGLAR (77)
Purity	98%
Source	Yeast
MW	9.1 kDa
Entrez	397072
Expiration Date	6 months from date of receipt.