

## Product Datasheet

### Ovine IL-4 Recombinant Protein (orb1820249)

<b>Description</b>	The Ovine IL-4 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The ...
<b>Target</b>	IL-4
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	HKCDITLLEEI IKTLNILTSR KNSCMELPVA DVFAAPKNAT EKETFCRAGI ELRRIYRSHM CLNKFLGGLD RNLSSLASKT CSVNEAKTST STLRDLLERL KTIMKEKYSK C (111)
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	12.6 kDa
<b>Entrez</b>	<a href="#">443327</a>
<b>Expiration Date</b>	6 months from date of receipt.