

Product Datasheet

Ovine IL-4 Recombinant Protein (orb1820250)

Description	The Ovine IL-4 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The ...
Target	IL-4
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	HKCDITLLEEI IKTLNILTSR KNSCMELPVA DVFAAPKNAT EKETFCRAGI ELRRIYRSHM CLNKFLGGLD RNLSSLASKT CSVNEAKTST STLRDLLERL KTIMKEKYSK C (111)
Purity	98%
Source	Yeast
MW	12.6 kDa
Entrez	443327
Expiration Date	6 months from date of receipt.