

Product Datasheet

Mouse ICAM-1 Recombinant Protein (orb1820252)

Description	The Mouse sICAM-1 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. T...
Target	ICAM-1
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	QVSIHPREAF LPQGGSVQVN CSSSCKEDLS LGLETQWLKD ELESGPNWKL FELSEIGEDS SPLCFENCGT VQSSASATIT VYSFPESVEL RPLPAWQQVG KDLTLRCHVD GGAPRTQLSA VLLRGEEILS RQPVGGHPKD PKEITFTVLA SRGDHGANFS CRTELDLRPQ GLALFSNVSE ARSLRTFDLP ATIPKLDTPD LLEVGTQQKL FCSLEGLFPA SEARIYLELG GQMPTQESTN SSDSVSATAL VEVTEFDRT LPLRCVLELA DQILETQRTL TVYNFSAPVL TLSQLEVSEG SQVTVKCEAH SGSKVLLSG VEPRPPTPQV QFTLNASSED HKRSFFCSAA LEVAGKFLFK NQTLELHVLY GPRLDETDCD GNWTWQEGSQ QTLKCAWGN PSPKMTERRK ADGALLPIGV VKSVKQEMNG TYVCHAFSSH GNVTRNVYLT VLYHSQNN (458)
Purity	98%
Source	Yeast
MW	50.2 kDa
Entrez	15894
Expiration Date	6 months from date of receipt.