

Product Datasheet

Bovine ICAM-1 Recombinant Protein (orb1820256)

Description	The Bovine sICAM-1 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. ...
Target	ICAM-1
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	GISIHPSKAI IPRGDSLTVN CSNSCDQKST FGLETVLIKE EVGRGDNWKV FQLRDVQEDI ELFCYSNCHK EQTIASMNLT VYWFPEHVEL APLPLWQPVG EELNLSCLVS GGAPRAHLSV VLLRGEEELG RQPVGKGEPA KVMFTVQSRR EDHGTNFSCR WELDLRSOGL ELFQNTSAPR KLQTYVLPSI DPHLEVPIIV EVGSRWPVNC TLDGLFPASD AKVYLVLGDQ KLESNITYDG DSVLAKAWME ENEEGTHSLK CSVTLGEVSR RTQENVTVYS FPLPTLTLSP PEVSEWTTVT VECVTRDGAV VKLNGTSAVP PGPRAQLKLN ASASDHRSNF SCSAALEIAG QVVKHKQTL LHVLYGPRLD QRDCPGNWTW QEGSEQTLKC EAQGNPIPKL NCSRKGDGAS LPIGDLRPVR REVAGTYLCR ATSARGRVTR EVVLNVLHGQ NIL (453)
Purity	98%
Source	Yeast
MW	50.0 kDa
Entrez	281839
Expiration Date	6 months from date of receipt.