

## Product Datasheet

### Bovine CX3CL1 Recombinant Protein (orb1820257)

<b>Description</b>	The Bovine CX3CL1 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. T...
<b>Target</b>	CX3CL1
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	QQHGVNKCNL FCNKMTQKIP ESRLVGYQRN RESCNDGAVI LKTVMGKSFV ADPKEEWVQK AMKHLDDHKAS VSQKSG (76)
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	8.6 kDa
<b>Entrez</b>	<b>517354</b>
<b>Expiration Date</b>	6 months from date of receipt.