

Product Datasheet

Bovine CX3CL1 Recombinant Protein (orb1820259)

Description	The Bovine CX3CL1 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. T...
Target	CX3CL1
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	QQHGVNKCNL FCNKMTQKIP ESRLVGYQRN RESCNDGAVI LKTVKGKSFC ADPKEEWVQK AMKHLDHKAS VSQKSG (76)
Purity	98%
Source	Yeast
MW	8.6 kDa
Entrez	517354
Expiration Date	6 months from date of receipt.