

Product Datasheet

Ferret IL-2 Recombinant Protein (orb1820263)

Description	The Ferret IL-2 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
Target	IL-2
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	APTTSSSTKE AQQQLEQLLL DLQLLLNGVK NYESPRMLTF KFYMPKKATE LTHLQCLAE LKLLEEVLYL AQSKNFHLTD IKELMSNINV TLLKLGSET SFKCEYDDET VTITEFLNKW ITFCQSIFST LT (132)
Purity	98%
Source	Yeast
MW	15.2 kDa
Entrez	101674110
Expiration Date	6 months from date of receipt.