

## Product Datasheet

# Canine TNFSF11 (RANKL) Recombinant Protein (orb1820267)

<b>Description</b>	The Canine TNFSF11 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. ...
<b>Target</b>	TNFSF11
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	EKAMMEGSWL EMARRGKTHT QPFAHLTINA TDIPSGSHKV SLSSWYHDRG WAKISNMTFS NGKLIVNQDG FYFLYANICF RHHETSGDLA TEYLQLMVYV TKTSIKIPSS HTLMKGGSTK YWSGNSEFHF YSINVGGFFK LRSGEEISIE VSNPSLLDPD QDATYFGAFK VLDID (175)
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	19.8 kDa
<b>Entrez</b>	<a href="#">609418</a>
<b>Expiration Date</b>	6 months from date of receipt.