

Product Datasheet

Canine TNFSF11 (RANKL) Recombinant Protein (orb1820269)

Description	The Canine TNFSF11 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. ...
Target	TNFSF11
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	EKAMMEGSWL EMARRGKTHT QPFAHLTINA TDIPSGSHKV SLSSWYHDRG WAKISNMTFS NGKLIVNQDG FYFLYANICF RHHETSGDLA TEYLQLMVYV TKTSIKIPSS HTLMKGGSTK YWSGNSEFHF YSINVGFFK LRSGEEISIE VSNPSLLDPD QDATYFGAFK VLDID (175)
Purity	98%
Source	Yeast
MW	19.8 kDa
Entrez	609418
Expiration Date	6 months from date of receipt.