

Product Datasheet

Canine IL-17F Recombinant Protein (orb1820272)

Description	The Canine IL-17F yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. T...
Target	IL-17F
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	ARKHLKAGET ALCPPLEDNS VRVDIRILRQ NRGISISNDF QNRSSSPWDY NITRDPHRFP SEIAEAQCRH SGCINAEGQE DSSMNSVPIQ QEFLVLRREP QGCSRSFRLE KVLVTVGCTC VTPIVRYVRA (130)
Purity	98%
Source	Yeast
MW	14.7 kDa
Entrez	481638
Expiration Date	6 months from date of receipt.