

Product Datasheet

Feline BAFF Recombinant Protein (orb1820273)

Description	The Feline BAFF yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
Target	BAFF
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	AALDPEDTVI QDCLQLIADS DTPTIRKGAY TFVPWLLSFK RGRALEEKEN KILVKETGYF FIYGQVLYTD NTFAMGHLIQ RKKVHVFGDE LSLVNLFRCI QNMPETLPNN SCYSAGIAKL EEGDELQLAI PREDAKISQD GDGTFFGALK LL (152)
Purity	98%
Source	Yeast
MW	17.1 kDa
Entrez	101088640
Expiration Date	6 months from date of receipt.