

Product Datasheet

Monkey IL-2 protein (orb1820277)

Description	The Monkey IL-2 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
Target	IL-2
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	APTSSSTKKT QLQLEHLLLD LQMILNGINN YKNPKLTRML TFKFYMPKKA TELRHLQCLE EELKPLEEVL NLAQSKSFHL RDTKDLISNI NVIVLELKGS ETTLMCEYAD ETATIVEFLN RWITFCQSII STLT (134)
Purity	98%
Source	Yeast
MW	15.5 kDa
Entrez	102129830
Expiration Date	6 months from date of receipt.