

## Product Datasheet

# Ferret IL-1 beta Recombinant Protein (orb1820278)

<b>Description</b>	The Ferret IL-1 beta yeast-derived recombinant protein is not tagged and is naturally endotoxin-free...
<b>Target</b>	IL-1 beta
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	APRSLDCKF QDINQKSLVL YNSYELRALH LNGTSVNQQA VFRMTFVQED EDITKIPVAL CIKEKNLYLS CVMKDGKPTL QLEMLDPKVY PKKRMEKRFV FNKTEIKKKV EFESSQFPNW YISTSQAEAM PVFLGNNRGG HDITDFTMEL SS (152)
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	17.6 kDa
<b>Entrez</b>	<a href="#">101688006</a>
<b>Expiration Date</b>	6 months from date of receipt.