

Product Datasheet

Ferret IL-1 beta Recombinant Protein (orb1820279)

Description	The Ferret IL-1 beta yeast-derived recombinant protein is not tagged and is naturally endotoxin-free...
Target	IL-1 beta
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	APRSLDCKF QDINQKSLVL YNSYELRALH LNGTSVNQQA VFRMTFVQED EDITKIPVAL CIKEKNLYLS CVMKDGKPTL QLEMLDPKVY PKKRMEKRFV FNKTEIKKKV EFESSQFPNW YISTSQAEAM PVFLGNNRGG HDITDFTMEL SS (152)
Purity	98%
Source	Yeast
MW	17.6 kDa
Entrez	101688006
Expiration Date	6 months from date of receipt.