

Product Datasheet

Dolphin IL-13 Recombinant Protein (orb1820283)

Description	The Dolphin IL-13 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. T...
Target	IL-13
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	GPVPPSTALK ELIEELVNIT QNQKAPLCNG SMVWSVNLTA NMYCSALEAL INVSNCSAIQ RTQRMLSALC LHKPSAGQVS SEHIRDTKIE VAHFLKDLLK HSRIIFHYGK FS (112)
Purity	98%
Source	Yeast
MW	12.4 kDa
Entrez	101338603
Expiration Date	6 months from date of receipt.