

Product Datasheet

Ferret IL-4 Recombinant Protein (orb1820291)

Description	The Ferret IL-4 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
Target	IL-4
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	HNFSIAIKEI IKTLNILTAR NDSCMELTVT EVFSAPKNTS EKEIFCRAAT VLQQLSTHNC SNRLLRGLHR NLRNMANMTC SVNEVKKSTL KDFLERLKVI MQQKYRHH (108)
Purity	98%
Source	Yeast
MW	12.5 kDa
Entrez	101686417
Expiration Date	6 months from date of receipt.