

Product Datasheet

Canine Fas Ligand Recombinant Protein (orb1820292)

Description	The Canine Fas Ligand yeast-derived recombinant protein is not tagged and is naturally endotoxin-free...
Target	Fas Ligand
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	QIQPNPPSE KRELKVAHL TGKPNRSIP LEWEDTYGIA LVSGVKYKKG GLVINDTGLY FVYSKVYFRG QSCNNKPLNH KVYMRNSKYP QDLMLMEGKI MNYCTTGQMW ARSSYLGAVF NLTSADHLYV NVSELSLVSF EESKTFFGLY KL (152)
Purity	98%
Source	Yeast
MW	17.3 kDa
Entrez	442968
Expiration Date	6 months from date of receipt.