

## Product Datasheet

# Canine Fas Ligand Recombinant Protein (orb1820294)

<b>Description</b>	The Canine Fas Ligand yeast-derived recombinant protein is not tagged and is naturally endotoxin-free...
<b>Target</b>	Fas Ligand
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	QIQPNPPSE KRELKVAHL TGKPNRSIP LEWEDTYGIA LVSGVKYKKG GLVINDTGLY FVYSKVYFRG QSCNNKPLNH KVYMRNSKYP QDLMLMEGKI MNYCTTGQMW ARSSYLGAVF NLTSADHLYV NVSELSLVSF EESKTFFGLY KL (152)
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	17.3 kDa
<b>Entrez</b>	<a href="#">442968</a>
<b>Expiration Date</b>	6 months from date of receipt.