

Product Datasheet

Guinea Pig IL-17A Recombinant Protein (orb1820307)

Description	The Guinea Pig IL-17A yeast-derived recombinant protein is not tagged and is naturally endotoxin-free...
Target	IL-17A
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	GIPIRNP GC PTATEGKNFL QNVKLNLSIF NPLTQNVNSR RSSDY YKRST SPWTLHRNEN PNRYPPIWE AECRYSGCVN AAGKEDHHVS SVPIQQEILV LQREPQNCPL SFRLEKMKVT VGCTCVTPIV RHVG (134) [Trp15Glu]
Purity	98%
Source	Yeast
MW	15.1 kDa
Entrez	100735572
Expiration Date	6 months from date of receipt.