

## Product Datasheet

# Guinea Pig IL-17A Recombinant Protein (orb1820308)

<b>Description</b>	The Guinea Pig IL-17A yeast-derived recombinant protein is not tagged and is naturally endotoxin-free...
<b>Target</b>	IL-17A
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	GIPIRNP GC PTATEGKNFL QNVKLNLSIF NPLTQNVNSR RSSDY YKRST SPWTLHRNEN PNRYPPIWE AECRYSGCVN AAGKEDHHVS SVPIQQEILV LQREPQNCPL SFRLEKMKVT VGCTCVPIV RHVG (134) [Trp15Glu]
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	15.1 kDa
<b>Entrez</b>	<a href="#">100735572</a>
<b>Expiration Date</b>	6 months from date of receipt.