

Product Datasheet

Caprine IFN gamma Recombinant Protein (orb1820318)

Description	The Caprine IFN gamma yeast-derived recombinant protein is not tagged and is naturally endotoxin-free...
Target	IFN gamma
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	QGPFFKEIEN LKEYFNASNP DVAKGGPLFS EILKNWKEES DKKIIQSQIV SFYFKLFENL KDNQVIQRSM DIIKQDMFQK FLNGSSEKLE DFKKLIQIPV DDLQIQRKAI NELIKVMNDL SPKSNLRKRK RSQNLFRGRR ASM (143)
Purity	98%
Source	Yeast
MW	16.9 kDa
Entrez	100860815
Expiration Date	6 months from date of receipt.