

Product Datasheet

Ferret IL-17A Recombinant Protein (orb1820319)

Description	The Ferret IL-17A yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. T...
Target	IL-17A
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	GIAFPQNPGC PNTEDKNFPQ HVKVNLNILN GNTNSRRPSD YHNRSTSPWN LHRNKDPERY PPVIWEAKCR HLGCVSADGT VNYHMNSVPI QQEILVLRRE SQHCPHSFRL EKMLVAVGCT CVTPIVRHVA (130)
Purity	98%
Source	Yeast
MW	14.8 kDa
Entrez	101672588
Expiration Date	6 months from date of receipt.