

## Product Datasheet

### Ferret IL-17A Recombinant Protein (orb1820321)

<b>Description</b>	The Ferret IL-17A yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. T...
<b>Target</b>	IL-17A
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	GIAFPQNPGC PNTEDKNFPQ HVKVNLNILN GNTNSRRPSD YHNRSTSPWN LHRNKDPERY PPVIWEAKCR HLGCVSADGT VNYHMNSVPI QQEILVLRRE SQHCPHSFRL EKMLVAVGCT CVTPIVRHVA (130)
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	14.8 kDa
<b>Entrez</b>	<b>101672588</b>
<b>Expiration Date</b>	6 months from date of receipt.