

Product Datasheet

Rat LIF Recombinant Protein (orb1820325)

Description	The Rat LIF yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The pur...
Target	LIF
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	SPLPITPVNA TCAIRHPCHG NLMNQIKSQL AQLNGSANAL FISYYTAQGE PFPNNVDKLC APNMTDFPPF HANGTEKTKL VELYRMVAYL GASLTNITWD QKLNPTAVS LQIKLNATTD VMRGLLSNVL CRLCNKYHVG HVDVPCVPDN SSKEAFQRKK LGCQLLGTYK QVISAVVQAF (180)
Purity	98%
Source	Yeast
MW	19.8 kDa
Entrez	60584
Expiration Date	6 months from date of receipt.