

## Product Datasheet

### Rat LIF Recombinant Protein (orb1820326)

|                         |  |
|-------------------------|--|
| <b>Description</b>      | The Rat LIF yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The pur...  |
| <b>Target</b>           | LIF  |
| <b>Form/Appearance</b>  | Lyophilized  |
| <b>Storage</b>          | -20°C  |
| <b>Note</b>             | For research use only  |
| <b>Protein Sequence</b> | SPLPITPVNA TCAIRHPCHG NLMNQIKSQL<br>AQLNGSANAL FISYYTAQGE PFPNNVDKLC<br>APNMTDFPPF HANGTEKTKL VELYRMVAYL<br>GASLTNITWD QKLNPTAVS LQIKLNATTD<br>VMRGLLSNVL CRLCNKYHVG HVDVPCVPDN<br>SSKEAFQRKK LGCQLLGTYK QVISAVVQAF<br>(180) |
| <b>Purity</b>           | 98%  |
| <b>Source</b>           | Yeast  |
| <b>MW</b>               | 19.8 kDa   |
| <b>Entrez</b>           | <b>60584</b>   |
| <b>Expiration Date</b>  | 6 months from date of receipt.   |