

## Product Datasheet

### Rat LIF Recombinant Protein (orb1820327)

<b>Description</b>	The Rat LIF yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The pur...
<b>Target</b>	LIF
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	SPLPITPVNA TCAIRHPCHG NLMNQIKSQL AQLNGSANAL FISYYTAQGE PFPNNVDKLC APNMTDFPPF HANGTEKTKL VELYRMVAYL GASLTNITWD QKLNPTAVS LQIKLNATTD VMRGLLSNVL CRLCNKYHVG HVDVPCVPDN SSKEAFQRKK LGCQLLGTYS QVISAVVQAF (180)
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	19.8 kDa
<b>Entrez</b>	<b>60584</b>
<b>Expiration Date</b>	6 months from date of receipt.