

Product Datasheet

Caprine LIF Recombinant Protein (orb1820330)

Description	The Caprine LIF yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
Target	LIF
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	SPLPITPVNA TCATRHPCPS NLMNQIRNQL GQLNSSANSL FILYYTAQGE PFPNNLDKLC SPNVTDFFPPF HANGTEKARL VELYRIIAYL GASLGNITRD QKVLNPYAHG LHSLKNATAD VLRGLLSNVL CRLCSKYHVS HVDVITYGPDY SGKDIFQKKK LGCQLLGKYK QVIQVLAQAF (180)
Purity	98%
Source	Yeast
MW	19.8 kDa
Entrez	102188524
Expiration Date	6 months from date of receipt.