

## Product Datasheet

# Cynomolgus Monkey LIF Recombinant Protein (orb1820331)

<b>Description</b>	The Cynomolgus Monkey Leukemia Inhibitory Factor (LIF) yeast-derived recombinant protein is not tagg...
<b>Target</b>	LIF
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	SPLPITPVNA TCAIRHPCHN NLMNQIRSQL AQLNGSANAL FILYYTAQGE PFPNNLDKLC GPNVTDFPPF HANGTEKAKL VELYRIVVYL GTSLGNITRD QKILNPSALS LHSKLNATAD ILRGLLSNVL CRLCSKYHVG HVDVTYGPDT SGKDFQKKK LGCQLLGKYK ETIAVLAQAF (180)
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	19.7 kDa
<b>Entrez</b>	<a href="#">102129162</a>
<b>Expiration Date</b>	6 months from date of receipt.