

Product Datasheet

Cynomolgus Monkey LIF Recombinant Protein (orb1820332)

Description	The Cynomolgus Monkey Leukemia Inhibitory Factor (LIF) yeast-derived recombinant protein is not tagg...
Target	LIF
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	SPLPITPVNA TCAIRHPCHN NLMNQIRSQL AQLNGSANAL FILYYTAQGE PFPNNLDKLC GPNVTDFPPF HANGTEKAKL VELYRIVVYL GTSLGNITRD QKILNPSALS LHSKLNATAD ILRGLLSNVL CRLCSKYHVG HVDVTYGPDT SGKDFVQKKK LGCQLLGKYK ETIAVLAQAF (180)
Purity	98%
Source	Yeast
MW	19.7 kDa
Entrez	102129162
Expiration Date	6 months from date of receipt.