

Product Datasheet

Dolphin LIF Recombinant Protein (orb1820335)

Description	The Dolphin Leukemia Inhibitory Factor (LIF) yeast-derived recombinant protein is not tagged and is ...
Target	LIF
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	SPLPITPVNA TCATRHPCHS NLMNQIKNQL AQLNSSAKAL FILYYTAQGE PFPNNLDKLC SPNVTDFFPPF HANGTERARL VELYRIIAYL GASLGNITRD QKVLNPSAHS LHSKLNATAD VMRGLLSNVL CRLCSKYHVA HVDVAYGPDY SGKDVVFQKKK LGCQLLGKYK QVIAVVAQAF (180) [Ser134Cys]
Purity	98%
Source	Yeast
MW	19.7 kDa
Entrez	101318472
Expiration Date	6 months from date of receipt.