

Product Datasheet

Mouse LIF Recombinant Protein (orb1820339)

Description	The Mouse Leukemia Inhibitory Factor (LIF) yeast-derived recombinant protein is not tagged and is na...
Target	LIF
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	SPLPITPVNA TCAIRHPCHG NLMNQIKNQL AQLNGSANAL FISYYTAQGE PFPNNVEKLC APNMTDFPSF HGNGTEKTKL VELYRMVAYL SASLTNITRD QKVLNPTAVS LQVKLNATID VMRGLLSNVL CRLCNKYRVG HVDVPPVPDH SDKEAFQRKK LGCQLLGTYK QVISVVVQAF (180)
Purity	98%
Source	Yeast
MW	19.9 kDa
Entrez	16878
Expiration Date	6 months from date of receipt.