

Product Datasheet

Caprine IL-6 Recombinant Protein (orb1820340)

Description	The Caprine IL-6 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. Th...
Target	IL-6
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	GPLGEDFKND TTPSRLLLT PEKTEALIKH IVDKISAIRK EICEKNDECE NSKETLAENK LKLPMEEKD GCFQSGFNQA ICLIKTTAGL LEYQIYLDLFL QNEFEGNQET VMELQSSIRT LIQILKEKIA GLITTPATNT DMLEKMQSSN EWWKNAKVII ILRSLENFLQ FSLRAIRMK (179)
Purity	98%
Source	Yeast
MW	20.4 kDa
Entrez	100860785
Expiration Date	6 months from date of receipt.