

Product Datasheet

Feline LIF Recombinant Protein (orb1820345)

Description	The Feline LIF yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The ...
Target	LIF
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	SPLPITPVNA TCATRHPCHS NLMNQIRNQL AQLNGSANAL FILYYTAQGE PFPNNLDKLC GPNVTDFPPF HANGTEKTRL VELYRIIAYL GASLGNITRD QKVLNPNALS LHSKLNATAD IMRGLLSNVL CRLCNKYHVA HVDVAYGPDT SGKDVVFQKKK LGCQLLGKYK QVIAVVAQAF [Met139Val](180)
Purity	98%
Source	Yeast
MW	19.7 kDa
Entrez	101090296
Expiration Date	6 months from date of receipt.