

Product Datasheet

Guinea Pig LIF Recombinant Protein (orb1820348)

Description	The Guinea Pig LIF yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. ...
Target	LIF
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	SPLPITPVNA TCATRHPCPS NVMSQIRNQL AQLNSSANAL FILYYTAQGE PFPNNLDKLC APNVTDFFPPF HANGTEKVKL VELYRIIAYL GAALGNITHD QKLLNPALS LHSKLSVTAD VMRGLLSNVH CRLCNKYHVA HVDVTYGPDT SDKDVFQKKK LGCQLLGKYK QVIGVVTQAF (180)
Purity	98%
Source	Yeast
MW	19.8 kDa
Entrez	100713602
Expiration Date	6 months from date of receipt.