

Product Datasheet

Swine Fas Ligand Recombinant Protein (orb1820349)

Description	The Swine Fas Ligand yeast-derived recombinant protein is not tagged and is naturally endotoxin-free...
Target	Fas Ligand
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	QIGHPNLPSE KKELRKVAHL TGKPNRSRIP LEWEDTYGIA LVSGVKYMKG SLVINDTGLY FVYSKVYFRG QYCNNQPLSH KVYTRNSRYP QDLVLMEGKM MNYCTTGQMW ARSSYLGAVF NLTSADHLYV NVSELSLVNF EESKTFGLY KL (152)
Purity	98%
Source	Yeast
MW	17.4 kDa
Entrez	396726
Expiration Date	6 months from date of receipt.