

Product Datasheet

Human LIF Recombinant Protein (orb1820354)

Description	The Human Leukemia Inhibitory Factor (LIF) yeast-derived recombinant protein is not tagged and is na...
Target	LIF
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	SPLPITPVNA TCAIRHPCHN NLMNQIRSQL AQLNGSANAL FILYYTAQGE PFPNNLDKLC GPNVTDFFPF HANGTEKAKL VELYRIVVYL GTSLGNITRD QKILNPSALS LHSKLNATAD ILRGLLSNVL CRLCSKYHVG HVDVTYGPDT SGKDVVFQKKK LGCQLLGKYK QIIAVLAQAF (180)
Purity	98%
Source	Yeast
MW	19.7 kDa
Entrez	3976
Expiration Date	6 months from date of receipt.