

Product Datasheet

Dolphin IL-2 Recombinant Protein (orb1820362)

Description	The Dolphin IL-2 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. Th...
Target	IL-2
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	APTSSSTENT KKQVQSLLQD LHLLLKEINN YENLKLFRML TFKFYMPKKA TELKHLQCLA EELKPLEDVL NVAQSKTQNS IDIKDLMDNI NRIVLTLKGS ETRFTCEYDD ETVTAVELLN KWITFCQSIY STMT (134)
Purity	98%
Source	Yeast
MW	15.5 kDa
Entrez	101325136
Expiration Date	6 months from date of receipt.