

## Product Datasheet

# Bovine Fas Ligand Recombinant Protein (orb1820367)

|                         |   |
|-------------------------|---|
| <b>Description</b>      | The Bovine Fas Ligand yeast-derived recombinant protein is not tagged and is naturally endotoxin-free...  |
| <b>Target</b>           | Fas Ligand  |
| <b>Form/Appearance</b>  | Lyophilized   |
| <b>Storage</b>          | -20°C   |
| <b>Note</b>             | For research use only   |
| <b>Protein Sequence</b> | QIGHPSPPSE KKEKKAHL TGKLNRSIP<br>LEWEDTYGIA LVSGVKYKKG<br>SLVINETGLY FVYSKVYFRG<br>QSCNNQPLSH KVYSRNFYYP<br>QDMVLMEGKM MNYCTSGKMW<br>ARSSYLGAVF NLTSADHLYV<br>NVSELSVSF EESKTFFGLY KL (152) |
| <b>Purity</b>           | 98%   |
| <b>Source</b>           | Yeast   |
| <b>MW</b>               | 17.3 kDa  |
| <b>Entrez</b>           | <a href="#">407111</a>  |
| <b>Expiration Date</b>  | 6 months from date of receipt.  |