

Product Datasheet

Bovine Fas Ligand Recombinant Protein (orb1820368)

Description	The Bovine Fas Ligand yeast-derived recombinant protein is not tagged and is naturally endotoxin-free...
Target	Fas Ligand
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	QIGHPSPPSE KKEKKAHL TGKLNRSIP LEWEDTYGIA LVSGVKYKKG SLVINETGLY FVYSKVYFRG QSCNNQPLSH KVYSRNFYYP QDMVLMEGKM MNYCTSGKMW ARSSYLGAVF NLTSADHLYV NVSELSLVSF EESKTFGLY KL (152)
Purity	98%
Source	Yeast
MW	17.3 kDa
Entrez	407111
Expiration Date	6 months from date of receipt.