

Product Datasheet

Rabbit CXCL3 Recombinant Protein (orb1820371)

Description	The Rabbit CXCL3 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. Th...
Target	CXCL3
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	APVVNELRCQ CLQTLQGIHL KNIQSVKVTT PGPHCDQTEV IASLKTGQEV CLNPTAPVVK KIIDKMLNKA SAN (73)
Purity	98%
Source	Yeast
MW	7.9 kDa
Entrez	100009192
Expiration Date	6 months from date of receipt.