

Product Datasheet

Canine CXCL13 Recombinant Protein (orb1820378)

Description	The Canine CXCL13 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. T...
Target	CXCL13
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	VLEAYNTNLK CKCIQEISGI IPIQLIERLQ IFPPGNGCPK REIIVWKKNK SVVCLNPIIK WTQKLNKILQ KKNATSTQPA PVLKKRIN (88)
Purity	98%
Source	Yeast
MW	10.0 kDa
Entrez	608156
Expiration Date	6 months from date of receipt.