

Product Datasheet

Caprine TNF alpha Recombinant Protein (orb1820386)

Description	The Caprine TNF alpha yeast-derived recombinant protein is not tagged and is naturally endotoxin-free...
Target	TNF alpha
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	LRSSSQASSN KPVAVHVANI SAPGQLRWGD SYANALKANG VELKDNQLVV PTDGLYLIYS QVLFGRHGCP STPLFLTHTI SRIAVSYQTK VNILSAIKSP CHRETPEGAE AKPWYEPIYQ GGVFQLEKGD RLSAEINQPE YLDYAESGQV YFGIIAL (157)
Purity	98%
Source	Yeast
MW	17.2 kDa
Entrez	100861232
Expiration Date	6 months from date of receipt.