

## Product Datasheet

# Caprine TNF alpha Recombinant Protein (orb1820388)

<b>Description</b>	The Caprine TNF alpha yeast-derived recombinant protein is not tagged and is naturally endotoxin-free...
<b>Target</b>	TNF alpha
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	LRSSSQASSN KPVAVVANI SAPGQLRWGD SYANALKANG VELKDNQLVV PTDGLYLIYS QVLFRGHGCP STPLFLTHTI SRIAVSYQTK VNILSAIKSP CHRETPEGAE AKPWYEPIYQ GGVFQLEKGD RLSAEINQPE YLDYAESGQV YFGIIAL (157)
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	17.2 kDa
<b>Entrez</b>	<a href="#">100861232</a>
<b>Expiration Date</b>	6 months from date of receipt.